




ORTHOPAEDICS INTERNATIONAL

Part of the Medical Aid International Family



External Fixation History and Training


History of External Fixation



**“he who desires
to practice
surgery must go
to war”**

Hippocrates c400BC


Ancient Egypt – 5000BC



Splints

Evidence found on
mummies
made from
bamboo, reeds and
wood or bark,
padded with linen

377BC - fracture splinting




Hippocrates first documented
the external fixator

1840 Jean-Francois Malgaigne




Credited with been first to produce an external fixation device using a crude pin

1893 - Keetley fixator



1894 - Clayton Parkhill




‘more
accurate
fixation of the
bones’

Bone Clamp

Wooden pins!


1902 - Freeman fixator



Belgian physician
Lambotte


Pins connected to each
other by an external
device - permitting
stabilization pins and
bone segments

1902 - Lambotte




First specially designed
device for external fixation

1934 - Roger Anderson




Developed an apparatus for the reduction of fractures which was made up of transfixion pins connected to metal clamps

1937 - Otto Strader



Developed a stabilization system for fractures that lead to a reduction of fractures in three planes


1938 - Raoul Hoffmann



Swiss Surgeon


First to patent
external fixation

1938 - Henry Jaquet



Forms a partnership with Raoul Hoffmann

1938 - Original design goals




3-dimensional

Universal ballpoint

Connected with rods


Reduction with fixator

1942 - Lewis and Breidenbach




First to describe the advantages of positioning the pins as far away as possible from the fracture
(near and far placement)

1958 - Vidal



Built on the Hoffman fixator designing a quadrilateral frame to guarantee fixation even more rigid and stable and demonstrating its utility in biomechanical studies


1964 - Gavril Ilizarov



Medical Aid
International

Supporting Healthcare in Low Resource Environments

1977 AO ASIF



Manual on AO/ASIF
Tubular External Fixator in
1977 released,
recommending the use of
external fixation. This
document provided
precise indications on use



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

1980 - DeBastiani



Verona

Orthofix

Monoaxial

Dynamic



Medical Aid
International

Supporting Healthcare in Low Resource Environments

1995 - Charles Taylor



Taylor Spatial frame

Computerised
distraction


Osteogenesis using an
octahedral assembly of
struts



Medical Aid
International

Supporting Healthcare in Low Resource Environments


1995 - Taylor Spacial Frame



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

External Fixation Training




**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

Four types of External Fixation

MONOLATERAL FIXATOR

- deformity correction
- limb lengthening



PIN-TO-BAR FIXATOR

- modular fixation in multiple planes
- temporary damage control



CIRCULAR FIXATOR

- for complex fractures and deformities in long bones




HYBRID FIXATOR

- for complex periarticular fractures in long bones



Types of fixator

- Monoplanar rigid frame
- Modular / Unilateral frames
- Bridging frames
- Peri-Articular frames
- Pin less Fixators
- Ring Fixators



For and against

Advantages

- ✓ Rapid
- ✓ Small incisions
- ✓ Easily removed
- ✓ Bone transport
- ✓ No knee pain
- ✓ 'Second hit'

Disadvantages

- ✗ Cumbersome
- ✗ Insufficient stability for weight bearing
- ✗ Loss of position
- ✗ Poor access to soft tissues
- ✗ Pin site infections
- ✗ Septic arthritis



Medical Aid
International

Supporting Healthcare in Low Resource Environments

Mechanisms of injury



www.UkF1.net



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

High energy impact



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments



Gunshot



Open Fracture



Soft Tissue



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

Why ex fix and not ORIF?

“First stage of staged treatment for challenging injuries”

- ✓ quick & minimally invasive fracture fixation
- ✓ allows soft-tissue recovery & stabilization
 - ✓ allows imaging of reduced fracture
- ✓ gives time for preparing a pre-operative plan for the final treatment
- ✓ voids complications & complex treatments



Indications

Impossible to nail

- ✓ too proximal or distal
- ✓ canal deformity
- ✓ skeletally immature
- ✓ battlefield trauma
- ✓ complex reconstruction

Rapid stabilisation

- ✓ soft tissues
- ✓ vascular injury
- ✓ peri-articular
- ✓ distraction device

Surgeon choice



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

Temporary spanning fixator



- ✓ Emergency stabilisation
- ✓ Soft Tissues are the Priority
- ✓ Poor soft tissues
- ✓ Vascular injury
- ✓ Compartment Syndrome



Medical Aid
International

Supporting Healthcare in Low Resource Environments

Diasphyseal/ Long bone fractures




Metaphyseal fractures



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments


Intra articular fractures



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments


Pelvic fractures



Medical Aid
International

Supporting Healthcare in Low Resource Environments

Wrist fracture



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

Bone transport



Distraction Osteogenesis

Bone in soft tissue sleeve

Large defects possible


Acute shortening with
later bone transport



Medical Aid
International

Supporting Healthcare in Low Resource Environments

Biomechanics



**Medical Aid
International**


Supporting Healthcare in Low Resource Environments

A.G.Apley

‘it joins two things together’




‘it needs a little movement’



Medical Aid
International

Supporting Healthcare in Low Resource Environments

Fracture healing - the basics



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

Principles - Pins



Conical
Cylindrical
Transfixion



Self Cutting/ Tapping
Blunt

Coatings - Hydroxyl Apatite
Silver Nitrite




Medical Aid
International

Supporting Healthcare in Low Resource Environments

Principles - Pins

Safe corridors

Pin and Wire placement



www.AONA.com



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments




Principles - Pins

Maximise **A**

Minimise **B**

Stack bars **C**

Additional planes **D**



Medical Aid
International

Supporting Healthcare in Low Resource Environments


Correct pin placement

A

At least 2 pins per segment

B

One pin allows bone movement around the pin



Medical Aid
International

Supporting Healthcare in Low Resource Environments

Correct pin placement

A


Spread the pins
'near and far'

B

Close together, pins
behave like one single
pin

B

A



Medical Aid
International

Supporting Healthcare in Low Resource Environments


Correct pin placement

A

Place frame close to the bone

B

A shorter lever arm reduces the momentum at the pin entry site





Medical Aid
International

Supporting Healthcare in Low Resource Environments

Pin placement

Penetrate the bone in the
middle

Best grip in the bone





Medical Aid
International

Supporting Healthcare in Low Resource Environments

Pin placement

- ✓ Travel the longest distance possible through the bone
- ✓ Bi-cortical purchase
- ✓ Avoid soft tissue damage as much as possible
- ✓ Think of patient comfort



Medical Aid
International

Supporting Healthcare in Low Resource Environments

Pin placement



- ✓ Pins should be 2 cm or more away from the fracture
- ✓ Insert the most difficult pins first
- ✓ Pay attention to where you can place pins:
 - joints
 - fracture
 - soft tissues



Medical Aid
International


Supporting Healthcare in Low Resource Environments

Elasticity and stability

A frame should be stable and
elastic

Stability avoids loss of reduction


Elasticity generates micro-
movements
thus callus



Medical Aid
International

Supporting Healthcare in Low Resource Environments


Frame stability - bookcase



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments


Frame stability - Ikea bookcase



Medical Aid
International

Supporting Healthcare in Low Resource Environments

Frame stability - Ikea bookcase




**Medical Aid
International**

Supporting Healthcare in Low Resource Environments


Pin placement

Extra rod increases
frame strength and
stability



Elasticity and pin diameter


Frame rigidity
(N/mm)



Medical Aid
International

Supporting Healthcare in Low Resource Environments

Rod stiffness



Medical Aid
International

Supporting Healthcare in Low Resource Environments


Pelvic stabilisation



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

Emergency stabilisation



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments

Thank you



 **Medical Aid International**
Supporting Healthcare in Low Resource Areas 

 **ORTHOPAEDICS
INTERNATIONAL**
Part of the Medical Aid International Family



**Medical Aid
International**

Supporting Healthcare in Low Resource Environments